

Supplementary Tables

Supplementary Table 1: Antibodies used for immunophenotyping in peripheral blood

PANEL	MARKER	CLONE	LABEL	ISOTYPE	PROVIDER
T cells	Dead/live	--	Violet	--	INVITROGEN
	CD3	UCHT1	FITC	IgG1	BECKMAN
	CD4	13B8.2	PC5.5	IgG1	BECKMAN
	CD8	B9.11	APCa750	IgG1	BECKMAN
	CD45RA	2H4	PC7	IgG1	BECKMAN
	CCR7	150503	PE	IgG2a	R&D Systems
	CD38	LS198-4-3	APC	IgG1	BECKMAN
	HLA-DR	Immu-357	KrO	IgG1	BECKMAN
Treg cells	Dead/live	--	Violet	--	INVITROGEN
	CD3	UCHT1	FITC	IgG1	BECKMAN
	CD4	13B8.2	PC5.5	IgG1	BECKMAN
	CD25	1HT44H3	APCa750	IgG2a	BECKMAN
	CD45RO	UCHL1	PC7	IgG1	BECKMAN
	CCR4	205410	PE	IgG2b	R&D Systems
	CD127	R34.34	APC	IgG1	BECKMAN
	HLA-DR	Immu-357	KrO	IgG1	BECKMAN
NK cells Monocytes Dendritic cells	Dead/live	--	Violet	--	INVITROGEN
	CD3	UCHT1	FITC	IgG1	BECKMAN
	CD19	J3-119	FITC	IgG1	BECKMAN
	CD20	B9E9 (HRC20)	FITC	IgG2a	BECKMAN
	CD123	SSDCLY107D2	PC5.5	IgG1	BECKMAN
	CD14	RMO52	APCa750	IgG2a	BECKMAN

	CD11c	BU15	PC7	IgG1	BECKMAN
	CD56	N901 (NKH.1)	PE	IgG1	BECKMAN
	CD16	3G8	APC	IgG1	BECKMAN
	HLA-DR	Immu-357	KrO	IgG1	BECKMAN